



<b>Equations &amp; Common Calculations</b>	<b>Identifier:</b> BPR - 0004
	<b>Revisions:</b> 0
	<b>Page:</b> 1 of 3

ISU Biological Sciences	Stockroom Procedure	Effective Date:03/19/2021
-------------------------	---------------------	---------------------------

## 1. Equations

1.1 **Molarity:** the amount of a substance in a certain volume of solution, moles of a solute per liters of a solution

1M = formula weight per liter (g/L)

Molar Mass = (g/mol)

$$M = \frac{\text{moles of solute}}{\text{liters of solution}}$$

**EX.** Prepare 600 mL of 0.4 M sucrose

Step 1. Find formula weight on the chemical bottle – F.W of sucrose is 342.3 g/mol

Step 2. Set up conversions, write out unit labels, and solve by cancelling units

$$\begin{array}{ccccccc} \text{Molar Mass} & & \times & \text{Molarity} & & \times & \text{Amount needed} \\ \\ \frac{342.3 \text{ g sucrose}}{1 \cancel{\text{ mol}}} & & \times & \frac{0.4 \cancel{\text{ mol}}}{1 \cancel{\text{ L}}} & & \times & 0.6 \cancel{\text{ L}} = 82.2 \text{ g sucrose} \end{array}$$

Add 82.2 grams of sucrose to 600 mL or 0.6 L of DH<sub>2</sub>O

1.2 **Dilutions:** decrease the concentration of a solute in a solution

$$C_1 \times V_1 = C_2 \times V_2 \quad \text{OR} \quad M_1 \times V_1 = M_2 \times V_2$$

**C<sub>1</sub>** is the initial (stock) concentration (molarity)      **V<sub>1</sub>** is the initial (stock) volume taken

**C<sub>2</sub>** is the concentration (molarity) of the dilution      **V<sub>2</sub>** is the final volume of the dilution



<b>Equations &amp; Common Calculations</b>	<b>Identifier: BPR - 0004</b> <b>Revisions: 0</b> <b>Page: 2 of 3</b>
--	---

ISU Biological Sciences	Stockroom Procedure	Effective Date: 03/19/2021
-------------------------	---------------------	----------------------------

1.3 **Percent Solutions:** amount or volume of chemical or compound per 100 mL of a solution

(Weight/Volume) = g per 100 mL of solvent

(Volume/volume) = mL per 100 mL of solvent

**EX.** Prepare a 5% NaHCO<sub>3</sub> (sodium bicarbonate) solution with a total volume of 500 mL

5% > 5g / 100 mL (w/v)

$$\frac{5 \text{ g NaHCO}_3}{100 \text{ mL DH}_2\text{O}} \times 500 \text{ mL} = 25 \text{ g NaHCO}_3$$

Add 25 g of sodium bicarbonate to 500 mL of DH<sub>2</sub>O



<b>Equations &amp; Common Calculations</b>	<b>Identifier:</b>	<b>BPR - 0004</b>
	<b>Revisions:</b>	<b>0</b>
	<b>Page:</b>	<b>3 of 3</b>

ISU Biological Sciences	Stockroom Procedure	Effective Date: 03/19/2021
-------------------------	---------------------	----------------------------

## 2. Common Calculations

### 2.1. CiDecon Disinfectant (1:128)

- 2.1.1. Prepare 3 gallons using the conversion 1 oz. concentrated CiDecon to 1 gallon DH<sub>2</sub>O
- 2.1.2. Fill gray spray bottles ¾ full and cap
- 2.1.3. Replace bi-weekly on Fridays in rooms 138, 139, 151, 152, 153, 243

### 2.2. Amphyl Disinfectant (1:200) or 1/2%

- 2.2.1. Prepare using the conversion table below
- 2.2.2. Pipette jars hold 6 liters (1.58 gallons) – add 31.6 mL Amphyl to 6 L DH<sub>2</sub>O
  - 2.2.2.1. 51 mL isopropyl alcohol/15 mL Amphyl can be added to maintain clarity if solution will not be used right away
  - 2.2.2.2. Use proper PPE – gloves, lab coat, and safety glasses
- 2.2.3. Refill once a semester

Dilution Strength	Amphyl	DH <sub>2</sub> O
(1:200)	5 mL	1 quart / 946.3 mL
(1:200)	20 mL	1 gallon / 3.785 L

### 2.3. Liquid Descaler / Acid Bath (1:124)

- 2.3.1. Prepare acid bath by using the conversion 8 mL descaler concentrate to 1 L DH<sub>2</sub>O
- 2.3.2. Use to soak glassware after decontamination, clean dishwasher, and clear rust
  - 2.3.2.1. Use proper PPE – gloves, lab coat, and safety glasses

### 2.4. Ethanol Dilutions from 190 proof / 95% ETOH

- 2.4.1. Use  $M_1 \times V_1 = M_2 \times V_2$ 
  - 2.4.1.1.  $M_1$  = the concentration of ethanol you want to prepare, ex. 70 %
  - 2.4.1.2.  $V_1$  = volume of ethanol you want to prepare
  - 2.4.1.3.  $M_2$  = the concentration of stock ethanol, 95 % or 190 proof
  - 2.4.1.4.  $V_2$  = volume of stock ethanol, 95 % or 190 proof
  - 2.4.1.5. Solve for  $V_2$ , subtract  $V_2$  from the total volume needed to determine final DH<sub>2</sub>O volume